


# Dr Christopher Orr Anterior Restoration Kit

Developed by Dr Christopher Orr to provide all the burs required for anterior restorations. This extensive kit includes the new Depth Marker Burs designed by Dr Zeynep Ayhan for use in predictable, precise and repeatable tooth preparation for porcelain laminate veneers, as well as all the burs needed for metal-ceramic and all-ceramic anterior crown preparations, refining anterior provisionals and finishing anterior ceramic restorations. Also included are burs for the placement and finishing of direct composites. The kit contains sixteen Schottlander Gold Diamond Burs consisting of one coarse, ten medium, four fine and one extra fine diamond burs.

<b>Initial Burs</b>						<b>Tooth Preparation</b>						<b>Other Applications</b>										
	<b>Note</b>	0.5mm	0.7mm	Round 1.0mm	Round 1.4mm		Round End Taper	<b>Note</b>	Round End Taper	Round End Cylinder	Round End Cylinder		Fine Flame	Rugby Ball	<b>Note</b>	Long Fine Flame	Short Fine Flame	Rugby Ball	Rugby Ball			
	<b>Use</b>	Depth Cuts	Depth Cuts	Depth Cuts	Depth Cuts		Remove old crown/veneer	<b>Use</b>	Veneer preparation	Crown prep axil reduction	Crown prep axil reduction		Break contacts	Palatal preparation	<b>Use</b>	Texture anterior composite	Refine margins on provisionals	Refine margins on provisionals	Fine adjustment of ceramic			
	<b>Quick Order Code</b>	M270	M271	M001	M002		X511	<b>Quick Order Code</b>	M320	M240	M244		M385	M052	<b>Quick Order Code</b>	F395	F397	F052	XF052			
<b>Refining the Preparation</b>						<b>Other Applications</b>					<b>Tooth Preparation</b>						<b>Refining the Preparation</b>					
	<b>Note</b>	Round End Taper	Rugby Ball	End Cutting	Round 1.0mm		Round 1.4mm	<b>Note</b>	Long Fine Flame	Short Fine Flame		Rugby Ball	Rugby Ball	<b>Note</b>	Round End Taper	Round End Cylinder		Round End Cylinder	Fine Flame	Rugby Ball		
	<b>Use</b>	Refine preparation	Refine preparation	Refine margins	Refine margins		Refine margins	<b>Use</b>	Texture anterior composite	Refine margins on provisionals		Refine margins on provisionals	Fine adjustment of ceramic	<b>Use</b>	Veneer preparation	Crown prep axil reduction		Crown prep axil reduction	Break contacts	Palatal preparation		
	<b>Quick Order Code</b>	F320	F052	M360	M001		M002	<b>Quick Order Code</b>	F395	F397		F052	XF052	<b>Quick Order Code</b>	M270	M271		M001	M002	X511		