

Dr Christopher Orr Posterior Restoration & Equilibration Kit

Developed by Dr Christopher Orr to provide the burs required for all posterior restorations (direct composite, all-ceramic inlays, onlays and crowns) as well as for occlusal equilibration. The kit contains fourteen burs: ten medium, two fine, one extra fine Schottlander Gold Diamonds and one Schottlander Fast-Cut CX Tungsten 2158 for efficient amalgam removal.

<p>Direct Restoration Preparation</p>			<p>Inlay/Onlay Preparation</p>											
	Note	Round End Cylinder		Tungsten CX Cross Cut	Note	Tungsten CX Cross Cut	Round End Taper	Long Fine Flame	Rugby Ball	Rugby Ball				
	Use	Caries access		Remove old amalgam	Use	Remove old amalgam	Internal preparation	Clear contacts	Cuspal reduction	Refine preparation				
	Quick Order Code	M296		2158	Quick Order Code	2158	M313	F395	M052	F052				
<p>Crown Preparation</p>									<p>Other Applications</p>					
	Note	Fine Flame	Round End Cylinder	Round End Cylinder	Round End Cylinder	Occlusal	Rugby Ball	Rugby Ball		End Cutting	Note	Equilibration Bur	Rugby Ball	Rugby Ball
	Use	Break contacts	Axial preparation	Axial preparation	Axial preparation	Occlusal surface preparation	Occlusal surface preparation	Refine preparation		Refine margins	Use	Adjust ceramic/occlusal equilibration	Adjust ceramic/provisionals	Fine adjustment of ceramic
	Quick Order Code	M385	M410	M240	M244	M041	M052	F052		M360	Quick Order Code	M262	F052	XF052