











Pattern Number	Schottlander Quick Order Code	Region of Use
Round End Taper 850 014 10 850 014 10	M320 F320	 Veneer and anterior crown preparation bur. Medium grit. Veneer and anterior crown preparation refining bur. Fine grit.
Round End Cylinder 881 010 8 881 010 8 (short shank) 881 014 8	M240 M241 M244	 Gold and ceramic preparation bur. Medium grit. Gold and ceramic preparation bur. Medium grit. Posterior PFM crown preparation bur. Medium grit.
Depth Cutter <i>(not included)</i> 834 0.5 (depth mm)	M270	 Veneer depth cutting bur. Medium grit.
End Cutting - 839 010	M360	 Refine margins. Medium grit.
Flame 859 010 10 863 014 10	M385 F395	 Break contacts. Medium grit. Anterior composite finishing bur. Fine grit.
Football 368 023 5 379 016 3.5 379 023 5	M052 F053 F054	 Occlusal adjustment. Medium grit. Fine occlusal adjustment bur. Fine grit. Denture adjustment bur. Occlusal rest seats. Fine grit.
Round - 801 018 1.8	F004	 Denture adjustment bur. Occlusal rest seats. Fine grit.
Flat End Cylinder - 835 012 5	M073	 Denture adjustment bur. Occlusal rest seats. Medium grit.
Cross Cut Domed Fissure - 012 4.5	GT42	 Cavity cutting.
Round End Taper 850 016 10	X511	 Remove crown/veneer.